

Rev05
Update: Nov,03,2023

DATASHEET

THE™ HA Tag Antibody [FITC], mAb, Mouse

Cat. No.: A01621

Overview

Specificity	GenScript THE™ HA Tag Antibody [FITC], mAb, Mouse recognizes C-terminal, N-terminal, and internal HA tagged fusion proteins.
Host Species	Mouse
Immunogen	A synthetic peptide YPYDVPDYA coupled to KLH
Conjugate	Fluorescein (FITC)

Applications

Working concentrations for specific applications should be determined by the investigator. The appropriate concentrations may be affected by secondary antibody affinity, antigen concentration, the sensitivity of the method of detection, temperature, the length of the incubations, and other factors. The suitability of this antibody for applications other than those listed below has not been determined. The following concentration ranges are recommended starting points for this product.


Application	Recommended Usage
Immunocytochemistry/Immunofluorescence (ICC/IF)	0.5-2 µg/ml
Flow Cytometry	0.5-2 µg for 1 x 10 ⁶ cells

Properties

Form	Lyophilized
Storage Buffer	lyophilized with PBS, pH 7.4, containing 1% BSA and 0.02% sodium azide
Reconstitution	Reconstitute the lyophilized powder with deionized water (or equivalent) to an final concentration of 1 mg/mL.
Storage Instructions	The lyophilized product remains stable for up to 1 year at -20 °C from the date of receipt. Upon reconstitution, it can be stored for 2-3 weeks at 2-8 °C or for up to 12 months at -20 °C or below. Avoid repeated freezing and thawing cycles.
Purification	Protein A affinity column

Isotype	Mouse IgG1, κ
Clonality	Monoclonal
Clone ID	5E11D8

Examples



Flow cytometric analysis of transfected HEK-293 cell expressing HA-tagged protein using THE™ Anti-HA-tag [FITC] Monoclonal Antibody (GenScript, A01621; shaded histogram) or with a negative control antibody (open histogram).

Background

Target Background : GenScript THE™ HA Tag Antibody [FITC], mAb, Mouse is FITC conjugated THE™ HA Tag Antibody, mAb, Mouse (A01244) purified from mice ascites by Protein A affinity column. This product is suitable for detecting the expression level of HA fusion proteins by various immunoassays. Absorption Wavelength: 492 nm Emission Wavelength: 518 nm Fluorochrome/Protein Ratio: >3 moles FITC per mole antibody

Synonyms : THE™ HA Tag Antibody [FITC], mAb, Mouse;

For laboratory research use only. Direct human use, including taking orally and injection and clinical use are forbidden.

生产商：南京金斯瑞生物科技有限公司 江苏省南京市江宁区科学园雍熙路28号

Manufacturer: Nanjing GenScript Biotech Co., Ltd. No. 28Yongxi Road, Jiangning District, Nanjing, Jiangsu, China