

DATASHEET

Version: 2016-08-17

GFP Antibody, pAb, Goat**Cat. No.:** A01703-40**Size:** 40 µg**Synonyms:** Green Fluorescent Protein antibody;**Description:**

Green fluorescence protein (GFP) is a 27 kDa protein derived from jellyfish *Aequorea victoria*. It emits green light (emission peak at a wavelength of 509 nm) when excited by blue light (excitation peak at a wavelength of 395 nm). GFP has become a very useful tool as a fusion protein to report gene expression, to trace cell lineage and to define subcellular protein localizations. YFP differs from GFP due to a mutation at T203Y.

GenScript GFP Antibody, pAb, Goat is developed in goat using purified recombinant full-length GFP protein.

Immunogen: purified recombinant full-length GFP protein**Host:** Goat**Conjugation:** Unconjugated**Formulation:**

0.5 mg/ml, lyophilized with PBS, pH 7.4, containing 0.02% sodium azide

Ig Subclass: IgG**Specificity:** **GenScript GFP Antibody, pAb, Goat**

specifically reacts with fusion proteins containing GFP epitope tag. This antibody recognizes all variants of GFP proteins such

as eGFP, AcGFP1 and cGFP. It slightly recognizes a variety of variants RFP proteins such as tdTomato, mOrange, DsRed, DsRed2 and mCherry.

Purification: Immunoaffinity chromatography**Applications:**

Working concentrations for specific applications should be determined by the investigator. The appropriate concentrations may be affected by secondary antibody affinity, antigen concentration, the sensitivity of the method of detection, temperature, the length of the incubations, and other factors. The suitability of this antibody for applications other than those listed below has not been determined. The following concentration ranges are recommended starting points for this product.

Western blot: 0.2-1 µg/ml**Immunoprecipitation:** 1-5 µg/ml**Other Applications:** user optimized**Reconstitution:**

Reconstitute the lyophilized powder with deionized water (or equivalent) to an antibody concentration of 0.5 mg/ml.

Storage:

The antibody is stable in lyophilized form if stored at -20°C or below. The reconstituted antibody can be stored for 2-3 weeks at 2-8°C. For long term storage, aliquot and store at -20°C or below. Avoid repeated freezing and thawing cycles.