

Rev09
Update: Sep,01,2026

DATASHEET

THE™ DYKDDDDK Tag Antibody [iFluor 488], mAb, Mouse

Cat. No.: A01809

Overview

Specificity	THE™ DYKDDDDK Tag Antibody [iFluor 488], mAb, Mouse recognizes N-terminal, internal and C-terminal Flag-tagged proteins.
Host Species	Mouse
Immunogen	A synthetic peptide DYKDDDDK coupled to KLH
Conjugate	iFluor 488

Applications

Working concentrations for specific applications should be determined by the investigator. The appropriate concentrations may be affected by antigen concentration, the sensitivity of the method of detection, temperature, the length of the incubations, and other factors. The suitability of this antibody for applications other than those listed below has not been determined. The following concentration ranges are recommended starting points for this product.

Application	Recommended Usage
Immunocytochemistry/Immunofluorescence (ICC/IF)	1-4 µg/ml
Flow Cytometry	1-4 µg/ml


Properties

Form	Lyophilized
Storage Buffer	lyophilized with PBS, pH 7.4, containing 10 mg/mL BSA and 0.02% sodium azide.
Reconstitution	Reconstitute the lyophilized powder with deionized water (or equivalent) to an final concentration of 0.5 mg/mL.
Storage Instructions	The lyophilized product remains stable for 2 years at -20°C from the date of receipt. Upon reconstitution, it can be stored for 2-3 weeks at 2-8°C or for 12 months at -20°C or below. Avoid repeated freezing and thawing cycles.
Purification	Protein A affinity column
Isotype	Mouse IgG2b, κ


GenScript USA, Inc.

Clonality	Monoclonal
Clone ID	5A8E5

Examples



Flow cytometric analysis of non-transfected CHO cells (Black) or Flag-tagged protein transfected CHO cells (Red) using **THE™ DYKDDDDK Tag Antibody [iFluor 488], mAb, Mouse** (A01809-100, 4 µg/ml).



Immunocytochemistry/Immunofluorescence analysis of non-transfected CHO cells (a) or Flag-tagged protein transfected CHO cells (b) using **THE™ DYKDDDDK Tag Antibody [iFluor 488], mAb, Mouse** (A01809-100, 4 µg/ml).

Background

Target Background : The Flag tag is a polypeptide of with the amino acid sequence DYKDDDDK amino acid sequence that can be added to a target protein using recombinant DNA technology. It can be fused to the N-terminus or C-terminus of the protein to facilitate detection and purification. An anti-Flag tag antibody is a useful tool for the analysis of Flag-tagged proteins. THE™ DYKDDDDK Tag Antibody [iFluor 488], mAb, Mouse is THE™ DYKDDDDK Tag Antibody (A00187) conjugated with iFluor 488 under optimal conditions with an F/P ratio of 4-8. It is suitable for detecting the expression level of DYKDDDDK-tagged proteins. iFluor 488 is a bright and photostable fluorescent dye. It is an excellent alternative to Alexa Fluor 488.

Synonyms : THE™ DYKDDDDK Tag Antibody [iFluor 488], mAb, Mouse; DYKDDDDK Tag Antibody [iFluor 488], mAb, Mouse; 488 conjugated Flag antibody; 488 conjugated Flag monoclonal antibody;

Terms of Use

Terms of Use**Third-Party Trademark Acknowledgement:**

All product and company names mentioned are trademarks or registered trademarks of their respective owners. iFluor® is a registered trademark of AAT Bioquest, Inc. and is used for identification purposes only. GenScript is not sponsored by, affiliated with, or endorsed by the respective trademark owners.

Limited Use License and Use Restrictions

Limited Use License and Use Restrictions

This Product is provided under a limited, non-exclusive, non-transferable, and non-sublicensable end-user license strictly limited to the permitted uses set forth in GenScript's official documentation, quotations, or written agreements. Any use beyond such scope requires GenScript's prior written authorization or license.

Except as expressly authorized in writing by GenScript, Customer shall not, and shall not permit, procure, enable, or assist any third party to: (a) commercialize or otherwise exploit the Product, including by reselling, distributing, transferring, sublicensing, or otherwise making the Product (including any product incorporating the Product) available to any third party; (b) copy, reproduce, clone, reverse clone, sequence, analyze, characterize, deconstruct, reconstruct, reverse engineer, or otherwise attempt to discover the composition, sequence, structure, binding characteristics, epitope recognition characteristics, or other proprietary features of the Product; (c) modify, derivatize, humanize, optimize, mutate, design around, or generate any variants, derivatives, fragments, substitutes, or functional equivalents of the Product; or (d) use the Product, any information derived therefrom, or any computational, AI-enabled, or sequence-based design method to develop, manufacture, offer for sale, sell, or commercialize any product or service competitive with the Product or GenScript's products or services.

All intellectual property rights and all rights not expressly granted herein are reserved by GenScript. No rights are granted by implication, estoppel, or otherwise.

For laboratory research use only. Direct human use, including taking orally and injection and clinical use are forbidden.

Manufacturer: Nanjing GenScript Biotech Co., Ltd. No. 28Yongxi Road, Jiangning District, Nanjing, Jiangsu, China