

Rev08
Update: Aug,31,2026

DATASHEET

Anti-Adalimumab Antibody, pAb, Rabbit

Cat. No.: A01886

Overview

Specificity	The product is specific for Adalimumab. This antibody serves as an excellent positive control for Adalimumab. The product is specific for Adalimumab. This antibody serves as an excellent positive control for Adalimumab immunogenicity (ADA) assays.
Host Species	Rabbit
Immunogen	Adalimumab
Conjugate	Unconjugated

Applications

Working concentrations for specific applications should be determined by the investigators. The appropriate concentration may be affected by secondary antibody affinity, antigen concentration, the sensitivity of the method of detection, temperature, the length of the incubations, and other factors. The suitability of this antibody for applications other than those listed below has not been determined. The following concentration ranges are recommended starting points for this product.

Application	Recommended Usage
ELISA	0.05-0.2 µg/ml

Properties

Form	Lyophilized
Storage Buffer	lyophilized with PBS, pH 7.4, containing 0.02% sodium azide.
Reconstitution	Reconstitute the lyophilized powder with deionized water (or equivalent) to an final concentration of 0.5 mg/mL.
Storage Instructions	The lyophilized product remains stable for 1 year at -20°C from the date of receipt. Upon reconstitution, it can be stored for 2-3 weeks at 2-8°C or for 12 months at -20°C or below. Avoid repeated freezing and thawing cycles.
Purification	Antigen affinity purification
Isotype	Rabbit IgG
Clonality	Polyclonal

GenScript USA, Inc.

Clone ID	Not applicable
Note	GenScript can customize this product per customer's request including product size, buffer components, etc.

Background

Target Background : Adalimumab (Humira) is a humanized monoclonal antibody that is approved by the U.S. Food and Drug Administration for the treatment of patients with rheumatoid arthritis, psoriatic arthritis, ankylosing spondylitis, Crohn's disease, ulcerative colitis, chronic psoriasis, hidradenitis suppurativa, and juvenile idiopathic arthritis. Adalimumab binds to the TNF α (tumor necrosis factor-alpha), and blocks its interaction with TNF α receptors, reducing the inflammatory response of autoimmune diseases. GenScript Anti-Adalimumab Antibody, pAb, Rabbit is produced from sera of a rabbit immunized with Adalimumab.

Synonyms : Rabbit polyclonal to Humira, Exemptia

Limited Use License and Use Restrictions

Limited Use License and Use Restrictions

This Product is provided under a limited, non-exclusive, non-transferable, and non-sublicensable end-user license strictly limited to the permitted uses set forth in GenScript' s official documentation, quotations, or written agreements. Any use beyond such scope requires GenScript' s prior written authorization or license.

Except as expressly authorized in writing by GenScript, Customer shall not, and shall not permit, procure, enable, or assist any third party to: (a) commercialize or otherwise exploit the Product, including by reselling, distributing, transferring, sublicensing, or otherwise making the Product (including any product incorporating the Product) available to any third party; (b) copy, reproduce, clone, reverse clone, sequence, analyze, characterize, deconstruct, reconstruct, reverse engineer, or otherwise attempt to discover the composition, sequence, structure, binding characteristics, epitope recognition characteristics, or other proprietary features of the Product; (c) modify, derivatize, humanize, optimize, mutate, design around, or generate any variants, derivatives, fragments, substitutes, or functional equivalents of the Product; or (d) use the Product, any information derived therefrom, or any computational, AI-enabled, or sequence-based design method to develop, manufacture, offer for sale, sell, or commercialize any product or service competitive with the Product or GenScript' s products or services.

All intellectual property rights and all rights not expressly granted herein are reserved by GenScript. No rights are granted by implication, estoppel, or otherwise.

For laboratory research use only. Direct human use, including taking orally and injection and clinical use are forbidden.

Manufacturer: Nanjing GenScript Biotech Co., Ltd. No. 28Yongxi Road, Jiangning District, Nanjing, Jiangsu, China