

Human IFN- γ Antibody (21G2), mAb, Mouse

Cat. No.: A01898

Overview

Specificity	The product is specific for human IFN- γ . This antibody binds with recombinant IFN- γ protein in ELISA and Western blot assay. The antibody is recommended as a capture antibody in a Sandwich ELISA assay with detection antibody GenScript, A01897-40, Human IFN- γ Antibody (16F3) [Biotin], mAb, Mouse.
Host Species	Mouse
Immunogen	Recombinant human IFN- γ
Conjugate	Unconjugated

Applications

Working concentrations for specific applications should be determined by the investigators. The appropriate concentration may be affected by secondary antibody affinity, antigen concentration, the sensitivity of the method of detection, temperature, the length of the incubations, and other factors. The suitability of this antibody for applications other than those listed below has not been determined. The following concentration ranges are recommended starting points for this product.


Application	Recommended Usage
ELISA	0.02-0.1 $\mu\text{g/ml}$
Western Blot	1 $\mu\text{g/ml}$

Properties


Form	Lyophilized
Storage Buffer	lyophilized with PBS, pH 7.4, containing 0.02% sodium azide.
Reconstitution	Reconstitute the lyophilized powder with deionized water (or equivalent) to an final concentration of 0.5 mg/mL.
Storage Instructions	The lyophilized product remains stable up to 1 year at -20 °C from date of receipt. Upon reconstitution, it can be stored for 2-3 weeks at 2-8 °C or for up to 12 months at -20 °C or below. Avoid repeated freezing and thawing cycles.

Purification	Protein A affinity column
Isotype	Mouse IgG1, κ
Clonality	Monoclonal
Clone ID	21G2
Note	GenScript can customize this product per customer's request including product size, buffer components, etc.

Examples



ELISA binding of Human IFN- γ Antibody (21G2), mAb, Mouse (GenScript, A01898-40) with Human IFN- γ recombinant protein, Coating antigen: Human IFN- γ recombinant protein, 1 μ g/ml, IFN- γ antibody dilution start from 1000 ng/ml, EC_{50} =15.36 ng/ml.



Standard curve of IFN- γ Sandwich ELISA. The IFN- γ Sandwich ELISA assay is developed by using Human IFN- γ Antibody(21G2), mAb, Mouse (GenScript, A01898-40) and Human IFN- γ Antibody (16F3) [Biotin], mAb, Mouse (GenScript, A01897-40) as the capture and detection antibodies, respectively.

The sensitivity of detecting IFN- γ is up to 1.56 ng/ml.

Background

Target Background : Interferon- γ (IFN- γ) is a pro-inflammatory cytokine that is central in host resistance to infection. It is mainly produced by natural killer cells and CD4+ and CD8+ T cells. Its receptors are found on nearly all cells, where it activates diverse responses that enable potential host cells to prevent invasive infection by bacteria, parasites, and viruses. Takayanagi et al. (2000) demonstrated that IFN- γ strongly suppresses osteoclastogenesis by interfering with the RANKL (602642)-RANK (603499) signaling pathway. Tsubota et al. (1999) reported that this upregulation in Sjogren syndrome patients may be controlled by interferon-gamma through the activation of transcription factor NF-KB.

Synonyms : Mouse monoclonal to IFNG,IFG,IFI,IFN-gamma,interferon gamma,Immune Interferon

For laboratory research use only. Direct human use, including taking orally and injection and clinical use are forbidden.