

Rev08
Update: Aug,31,2026

DATASHEET

Anti-Rituximab Antibody (17B6)[Biotin], mAb, Mouse

Cat. No.: A01971

Overview

Specificity	The product is specific for Rituximab. The antibody is recommended as a detection antibody in a pharmacokinetic (PK) bridging assay with capture antibody GenScript, A01969-40, Anti-Rituximab Antibody (6C1), mAb, Mouse.
Host Species	Mouse
Immunogen	Rituximab
Conjugate	Biotin

Applications

Working concentrations for specific applications should be determined by the investigators. The appropriate concentrations may be affected by secondary antibody affinity, antigen concentration, the sensitivity of the method of detection, temperature, the length of the incubations, and other factors. The suitability of this antibody for applications other than those listed below has not been determined. The following concentration ranges are recommended starting points for this product

Application	Recommended Usage
ELISA	0.01-0.1 µg/ml
Sandwich ELISA	A01969 (C)-A01971 (D)


Properties

Form	Lyophilized
Storage Buffer	lyophilized with PBS, pH 7.4, containing 1% BSA and 0.02% sodium azide.
Reconstitution	Reconstitute the lyophilized powder with deionized water (or equivalent) to an final concentration of 0.5 mg/mL.
Storage Instructions	The lyophilized product remains stable for 1 year at -20°C from the date of receipt. Upon reconstitution, it can be stored for 2-3 weeks at 2-8°C or for 12 months at -20°C or below. Avoid freeze/thaw cycles.

Purification	Protein A affinity column
Isotype	Mouse IgG1
Clonality	Monoclonal
Clone ID	17B6
Note	GenScript can customize this product per customer's request including product size, buffer components, etc.


Examples

Biotin conjugated Anti-Rituximab Antibody (17B6) binds with Rituximab



ELISA binding of Anti-Rituximab Antibody (17B6) [Biotin] (GenScript, A01971-40) with Rituximab. While the antibody does not recognize the human IgG (data not shown). Coating antigen: Rituximab, 1 µg/ml. Anti-Rituximab Antibody (17B6) (GenScript, A01971-40) dilutions start from 2,000 ng/ml. $EC_{50} = 5.075$ ng/ml.

Standard curve of Rituximab Sandwich ELISA



Standard curve of Rituximab Sandwich ELISA. The Rituximab Sandwich ELISA assay is developed by using Anti-Rituximab Antibody (6C1), mAb, Mouse (GenScript, A01969-40) and Anti-Rituximab Antibody (17B6), mAb, Mouse (GenScript, A01971-40) as the capture and detection antibodies, respectively. The sensitivity is less than 1 ng/ml.

Background

Target Background : Rituximab (Rituxan) is a humanized monoclonal antibody that is approved by the U.S. Food and Drug Administration for the treatment of hematological cancers and autoimmune diseases. When Rituximab binds to CD20 on the surface of normal and malignant B cells, natural killer (NK) cells is activated to destroy these B cells. It is therefore used to treat diseases with too many B cells, overactive B cells, or dysfunctional B cells. Anti-Rituximab Antibody (17B6), mAb, Mouse is produced from a hybridoma resulting from the fusion of partner and B-lymphocytes obtained from a mouse immunized with Rituximab.

Synonyms : Mouse monoclonal to Rituxan

Limited Use License and Use Restrictions

Limited Use License and Use Restrictions

This Product is provided under a limited, non-exclusive, non-transferable, and non-sublicensable end-user license strictly limited to the permitted uses set forth in GenScript's official documentation, quotations, or written agreements. Any use beyond such scope requires GenScript's prior written authorization or license.

Except as expressly authorized in writing by GenScript, Customer shall not, and shall not permit, procure, enable, or assist any third party to: (a) commercialize or otherwise exploit the Product, including by reselling, distributing, transferring, sublicensing, or otherwise making the Product (including any product incorporating the Product) available to any third party; (b) copy, reproduce, clone, reverse clone, sequence, analyze, characterize, deconstruct, reconstruct, reverse engineer, or otherwise attempt to discover the composition, sequence, structure, binding characteristics, epitope recognition characteristics, or other proprietary features of the Product; (c) modify, derivatize, humanize, optimize, mutate, design around, or generate any variants, derivatives, fragments, substitutes, or functional equivalents of the Product; or (d) use the Product, any information derived therefrom, or any computational, AI-enabled, or sequence-based design method to develop, manufacture, offer for sale, sell, or commercialize any product or service competitive with the Product or GenScript's products or services.

All intellectual property rights and all rights not expressly granted herein are reserved by GenScript. No rights are granted by implication, estoppel, or otherwise.

For laboratory research use only. Direct human use, including taking orally and injection and clinical use are forbidden.

Manufacturer: Nanjing GenScript Biotech Co., Ltd. No. 28Yongxi Road, Jiangning District, Nanjing, Jiangsu, China