

Rev09
Update: Aug,31,2026

DATASHEET

MonoRab™ Anti-Nivolumab Antibody (40F11) [Biotin], mAb, Rabbit

Cat. No.: A01989

Overview

Specificity	The antibody is recommended as a detection antibody in a pharmacokinetic (PK) bridging assay with capture antibody GenScript, A01987-40, MonoRab Anti-Nivolumab Antibody (33H4), mAb, Rabbit.
Host Species	Rabbit
Immunogen	Nivolumab
Conjugate	Biotin

Applications

Working concentrations for specific applications should be determined by the investigators. The appropriate concentrations may be affected by secondary antibody affinity, antigen concentration, the sensitivity of the method of detection, temperature, the length of the incubations, and other factors. The suitability of this antibody for applications other than those listed below has not been determined. The following concentration ranges are recommended starting points for this product.

Application	Recommended Usage
Sandwich ELISA	Capture: 0.01-0.1 µg/ml
Competitive ELISA	10-50 µg/ml


Properties

Form	Lyophilized
Storage Buffer	Lyophilized with PBS, pH 7.4, containing 1% BSA and 0.02% sodium azide.
Reconstitution	Reconstitute the lyophilized powder with deionized water (or equivalent) to a final concentration of 0.5 mg/ml.
Storage Instructions	The lyophilized product remains stable for 1 year at -20°C from the date of receipt. Upon reconstitution, it can be stored for 2-3 weeks at 2-8°C or for 12 months at -20°C or below. Avoid repeated freezing and thawing cycles.
Purification	Protein A affinity column

GenScript USA, Inc.

Isotype	Rabbit IgG
Clonality	Monoclonal
Clone ID	40F11
Note	GenScript can customize this product per customer's request including product size, buffer components, etc.

Examples




ELISA binding of MonoRab™ Anti-Nivolumab Antibody (40F11) [Biotin] (GenScript, A01989-40) with Nivolumab. While the antibody does not recognize the human IgG Fc fragment (data not shown).

Coating antigen: Nivolumab, 1 µg/ml.

MonoRab™ Anti-Nivolumab Antibody (40F11)[Biotin] (GenScript, A01989-40) dilutions start from 1,000 ng/ml.

EC₅₀= 14.30 ng/ml.



Standard curve of Nivolumab Sandwich ELISA. The Nivolumab Sandwich ELISA assay is developed by using MonoRab™ Anti-Nivolumab Antibody (33H4), mAb, Rabbit (GenScript, A01987-40) and MonoRab™ Anti-Nivolumab Antibody (40F11)[Biotin] (Genscript, A01989-40) as the capture and detection antibodies, respectively.

The sensitivity of detecting Nivolumab is up to 78 pg/ml.

Background

Target Background : Nivolumab is a fully human monoclonal IgG4 antibody that is approved by the U.S. Food and Drug Administration for the treatment of bladder cancer. It binds to the programmed cell death protein 1 (PD-1) on the surface of activated T cells. Nivolumab functions as PD-1 inhibitor for targeted immunotherapy. GenScript MonoRab™ Anti-Nivolumab Antibody (40F11), mAb, Rabbit is produced from a hybridoma resulting from the fusion of partner and B-lymphocytes obtained from a rabbit immunized with Nivolumab.

Synonyms : Rabbit monoclonal to Opdivo

Limited Use License and Use Restrictions

Limited Use License and Use Restrictions

This Product is provided under a limited, non-exclusive, non-transferable, and non-sublicensable end-user license strictly limited to the permitted uses set forth in GenScript's official documentation, quotations, or written agreements. Any use beyond such scope requires GenScript's prior written authorization or license.

Except as expressly authorized in writing by GenScript, Customer shall not, and shall not permit, procure, enable, or assist any third party to: (a) commercialize or otherwise exploit the Product, including by reselling, distributing, transferring, sublicensing, or otherwise making the Product (including any product incorporating the Product) available to any third party; (b) copy, reproduce, clone, reverse clone, sequence, analyze, characterize, deconstruct, reconstruct, reverse engineer, or otherwise attempt to discover the composition, sequence, structure, binding characteristics, epitope recognition characteristics, or other proprietary features of the Product; (c) modify, derivatize, humanize, optimize, mutate, design around, or generate any variants, derivatives, fragments, substitutes, or functional equivalents of the Product; or (d) use the Product, any information derived therefrom, or any computational, AI-enabled, or sequence-based design method to develop, manufacture, offer for sale, sell, or commercialize any product or service competitive with the Product or GenScript's products or services.

All intellectual property rights and all rights not expressly granted herein are reserved by GenScript. No rights are granted by implication, estoppel, or otherwise.

For laboratory research use only. Direct human use, including taking orally and injection and clinical use are forbidden.

Manufacturer: Nanjing GenScript Biotech Co., Ltd. No. 28Yongxi Road, Jiangning District, Nanjing, Jiangsu, China