

Rev06
Update: Aug,11,2025

DATASHEET

MonoRab™ Rabbit Anti-Camelid VHH Cocktail [PE]

Cat. No.: A02018

Overview

Specificity	GenScript MonoRab™ Rabbit Anti-Camelid VHH Cocktail [PE] is specific for Camelid VHH antibody (Llama, Alpaca and Camel). It has no cross-reactivity with mouse, rat, rabbit, goat or human immunoglobulins.
Host Species	Rabbit
Immunogen	Llama VHH antibody
Conjugate	Phycoerythrin (PE)

Applications

Working concentrations for specific applications should be determined by the investigators. The appropriate concentrations may be affected by primary/secondary antibody affinity, antigen concentration, the sensitivity of the method of detection, temperature, the length of the incubations, and other factors. The suitability of this antibody for applications other than those listed below has not been determined. The following concentration ranges are recommended starting points for this product. Optimal working dilutions must be determined by the end user.


Application	Recommended Usage
Flow Cytometry	1:50-1:800

Properties

Form	Liquid
Storage Buffer	Supplied in PBS (pH 7.4), containing gelatin and 0.03% ProClin300.
Concentration	0.5 mg/ml
Storage Instructions	Stored at 4°C and it should be protected from prolonged exposure to light. Do not freeze. This product is stable for 1 year upon receipt, when handled and stored as Instructed.
Purification	Protein A affinity column
Isotype	Rabbit IgG
Clonality	Cocktail

Note GenScript can customize this product per customer's request including product size, buffer components, etc.

Examples



FACS binding of a series of dilution of MonoRab™ Rabbit Anti-Camelid VHH Cocktail [PE] (Cat.no A02018) with Jurkat cells and VHH-based Jurkat-CAR cells.

Background

Target Background : A camelid single-domain antibody (sdAb, called Nanobody by Ablynx, the developer) is a peptide chain of about 110 amino acids long, comprising one variable domain of a heavy-chain antibody, these are called VHH fragments. Like a whole antibody, it is able to bind selectively to a specific antigen. Single-domain antibodies allow a broad range of applications in biotechnical as well as therapeutic uses due to their small size, simple production and high affinity. GenScript MonoRab™ Rabbit Anti-Camelid VHH Cocktail [PE] is PE conjugated MonoRab™ Rabbit Anti-Camelid VHH Cocktail. The Cocktail is a unique secondary antibody, composed of several rabbit anti-Camelid VHH monoclonal antibodies.

Synonyms : PE conjugated rabbit anti-camelid VHH, single-domain antibody monoclonal antibody cocktail (Min X Ms, Rt, Rb, Gt, Hu)

For laboratory research use only. Direct human use, including taking orally and injection and clinical use are forbidden.

Manufacturer: Nanjing GenScript Biotech Co., Ltd. No. 28Yongxi Road, Jiangning District, Nanjing, Jiangsu, China