

Rev09
Update: Aug,31,2026

DATASHEET

Anti-Denosumab Antibody (9A10), mAb, Mouse

Cat. No.: A02102

Overview

Specificity	The product is specific for Denosumab
Host Species	Mouse
Immunogen	Denosumab
Conjugate	Unconjugated

Applications

Working concentrations for specific applications should be determined by the investigators. The appropriate concentrations may be affected by secondary antibody affinity, antigen concentration, the sensitivity of the method of detection, temperature, the length of the incubations, and other factors. The suitability of this antibody for applications other than those listed below has not been determined. The following concentration ranges are recommended starting points for this product.

Application	Recommended Usage
ELISA	0.001-0.05 µg/ml
Competitive ELISA	5-30 µg/ml


Properties

Form	Lyophilized
Storage Buffer	This product is lyophilized with PBS, pH 7.2, containing 0.02% sodium azide.
Reconstitution	Reconstitute the lyophilized powder with deionized water (or equivalent) to a final concentration of 0.5 mg/ml.
Storage Instructions	The lyophilized product remains stable for 1 year at -20°C from the date of receipt. Upon reconstitution, it can be stored for 2-3 weeks at 2-8°C or for 12 months at -20°C or below. Avoid repeated freezing and thawing cycles.
Purification	Protein A affinity column
Isotype	Mouse IgG1, κ
Clonality	Monoclonal

Clone ID	9A10
Note	GenScript can customize this product per customer's request including product size, buffer components, etc.

Examples

Anti-Denosumab Antibody (9A10) blocks Denosumab binding with Human RANK ligand



Anti-Denosumab Antibody (9A10) (GenScript, A02102-40) blocks Denosumab binding with Human RANKL recombinant protein.


Coating antigen: Denosumab 1 µg/ml.

RANKL final concentration: 12.5 ng/ml.

Anti-Denosumab antibody dilutions start from 30 µg/ml.

IC₅₀ = 0.893 µg/ml.

Anti-Denosumab Antibody (9A10) binds with Denosumab




ELISA binding of Anti-Denosumab Antibody (9A10), mAb, Mouse (GenScript, A02102-40) with Denosumab. While the antibody does not recognize the human IgG (data not shown).

Coating antigen: Denosumab, 1 µg/ml.

Anti-Denosumab Antibody (9A10), mAb, Mouse (GenScript, A02102-40) dilution start from 1,000 ng/ml.

EC₅₀ = 2.921 ng/ml.

Standard curve of Denosumab Sandwich ELISA



Standard curve of Denosumab Sandwich ELISA. The Denosumab Sandwich ELISA assay is developed by using Anti-Denosumab Antibody (2H7), mAb, Mouse (GenScript, A02101-40) and Anti-Denosumab Antibody (9A10), mAb, Mouse (GenScript, A02102-40) as the capture and detection antibodies, respectively.

In this ELISA assay, Anti-Denosumab Antibody (9A10), mAb, Mouse (GenScript, A02102-40) was labeled with Biotin.

GenScript can provide customized conjugation services for this product per the customer's request.

The sensitivity of detecting Denosumab is up to 0.26 ng/ml.

Background

Target Background : Denosumab with trade name Prolia, is an FDA-approved drug for the treatment of patients with osteoporosis, treatment-induced bone loss, metastases to bone, and giant cell tumor of bone. Denosumab is a fully human monoclonal antibody that works as a RANK ligand (RANKL) inhibitor. Denosumab can bind RANKL, preventing RANKL from activating RANK, its receptor on the osteoclast surface. Anti-Denosumab Antibody (9A10), mAb, Mouse is produced from a hybridoma resulting from the fusion of partner and B-lymphocytes obtained from a mouse immunized with Denosumab.

Synonyms : Mouse monoclonal to Denosumab/Prolia

Limited Use License and Use Restrictions

Limited Use License and Use Restrictions

This Product is provided under a limited, non-exclusive, non-transferable, and non-sublicensable end-user license strictly limited to the permitted uses set forth in GenScript's official documentation, quotations, or written agreements. Any use beyond such scope requires GenScript's prior written authorization or license. Except as expressly authorized in writing by GenScript, Customer shall not, and shall not permit, procure, enable, or assist any third party to: (a) commercialize or otherwise exploit the Product, including by reselling, distributing, transferring, sublicensing, or otherwise making the Product (including any product incorporating the Product) available to any third party; (b) copy, reproduce, clone, reverse clone, sequence, analyze, characterize, deconstruct, reconstruct, reverse engineer, or otherwise attempt to discover the composition, sequence, structure, binding characteristics, epitope recognition characteristics, or other proprietary features of the Product; (c) modify, derivatize, humanize, optimize, mutate, design around, or generate any variants, derivatives, fragments, substitutes, or functional equivalents of the Product; or (d) use the Product, any information derived therefrom, or any computational, AI-enabled, or sequence-based design method to develop, manufacture, offer for sale, sell, or commercialize any product or service competitive with the Product or GenScript's products or services. All intellectual property rights and all rights not expressly granted herein are reserved by GenScript. No rights are granted by implication, estoppel, or otherwise.

For laboratory research use only. Direct human use, including taking orally and injection and clinical use are forbidden.

Manufacturer: Nanjing GenScript Biotech Co., Ltd. No. 28Yongxi Road, Jiangning District, Nanjing, Jiangsu, China