

Rev05
Update: Jan,07,2025

DATASHEET

Anti-Ustekinumab Antibody, pAb, Rabbit

Cat. No.: A02123

Overview

Specificity	This product is specific for Ustekinumab
Host Species	Rabbit
Immunogen	Ustekinumab
Conjugate	Unconjugated

Applications

Working concentrations for specific applications should be determined by the investigators. The appropriate concentrations may be affected by secondary antibody affinity, antigen concentration, the sensitivity of the method of detection, temperature, the length of the incubations, and other factors. The suitability of this antibody for applications other than those listed below has not been determined. The following concentration ranges are recommended starting points for this product.

Application	Recommended Usage
ELISA	0.001-0.1 µg/ml


Properties

Form	Lyophilized
Storage Buffer	This product is lyophilized with PBS, pH 7.2, containing 0.02% sodium azide.
Reconstitution	Reconstitute the lyophilized powder with deionized water (or equivalent) to a final concentration of 0.5 mg/ml.
Storage Instructions	The lyophilized product remains stable up to 1 year at -20 °C from date of receipt. Upon reconstitution, it can be stored for 2-3 weeks at 2-8 °C or for up to 12 months at -20 °C or below. Avoid repeated freezing and thawing cycles.
Purification	Antigen affinity purification
Isotype	Rabbit IgG
Clonality	Polyclonal

Note GenScript can customize this product per customer's request including product size, buffer components, etc.

Examples

Anti-Ustekinumab Antibody, pAb, Rabbit binds with Ustekinumab



ELISA binding of Anti-Ustekinumab Antibody, pAb, Rabbit (GenScript, A02123-40) with Ustekinumab.

Coating antigen: Ustekinumab, 1 µg/ml.

Anti-Ustekinumab Antibody, pAb, Rabbit (GenScript, A02123-40) dilution start from 1,000 ng/ml.

EC₅₀= 9.187 ng/ml.

Background

Target Background : Ustekinumab with trade name Stelara, is an FDA-approved drug to Crohn's disease, Ulcerative Colitis, Plaque Psoriasis and Psoriatic Arthritis, targeting both interleukins (IL)-12 and IL-23. Ustekinumab is a fully human immunoglobulin G1 kappa mAb that binds to IL-12 and IL-23. Ustekinumab binds to the p40 subunit common to IL-12 and IL-23 and prevents their interaction with the IL-12 receptor β1 subunit of the IL-12 and IL-23 receptor complexes, thereby disrupting the proinflammatory pathway. GenScript Anti-Ustekinumab Antibody, pAb, Rabbit is produced from sera of a rabbit immunized with Ustekinumab.

Synonyms : Rabbit polyclonal to Ustekinumab/Stelara

For laboratory research use only. Direct human use, including taking orally and injection and clinical use are forbidden.

Manufacturer: Nanjing GenScript Biotech Co., Ltd. No. 28Yongxi Road, Jiangning District, Nanjing, Jiangsu, China