

Rev01
Update: May,17,2023

DATASHEET

MonoRab™ Rabbit Anti-Humanized VHH Antibody [PE], mAb

Cat. No.: A02171

Overview

Specificity	The product is specific for Humanized VHH antibodies and Camelid (Llama, Alpaca, and Camel) heavy chain antibodies (IgG2&3), but does not recognize the camelid conventional antibody (IgG1).
Host Species	Rabbit
Immunogen	VHH antibody
Species Reactivity	It has no cross-reactivity with mouse, rat, rabbit, or goat immunoglobulins.
Conjugate	PE

Applications

Working concentrations for specific applications should be determined by the investigators. The appropriate concentrations may be affected by secondary antibody affinity, antigen concentration, the sensitivity of the method of detection, temperature, the length of the incubations, and other factors. The suitability of this antibody for applications other than those listed below has not been determined. The following concentration ranges are recommended starting points for this product.

Application	Recommended Usage
Flow Cytometry	1:100-1:1000

Properties

Form	Liquid
Storage Buffer	Supplied in PBS (pH 7.2), containing gelatin and 0.03% ProClin300.
Concentration	0.5 mg/ml
Storage Instructions	Stored at 4°C and it should be protected from prolonged exposure to light. Do not freeze. This product is stable for 1 year upon receipt, when handled and stored as Instructed.
Purification	Protein A affinity column

Isotype	Rabbit IgG
Clonality	Monoclonal
Clone ID	45H8
Note	GenScript can customize this product per customer's request including product size, buffer components, etc.

Background

Target Background : A camelid single-domain antibody (sdAb, called Nanobody by Ablynx, the developer) is a peptide chain of about 110 amino acids long, comprising one variable domain of a heavy-chain antibody, these are called VHH fragments. Like a whole antibody, it is able to bind selectively to a specific antigen. Single-domain antibodies allow a broad range of applications in biotechnical as well as therapeutic uses due to their small size, simple production, and high affinity. However, their sequences may increase the risk of immunogenicity and anti-drug antibody (ADA) development in humans, and thus, sdAbs are routinely humanized during development.

Synonyms : Single-domain antibody, sdAb, Variable Domain of Heavy Chain of Heavy-Chain Antibody, Nanobody

Reference : 1. [Immunogenicity and humanization of single-domain antibodies.](#)
2. [Properties, production, and applications of camelid single-domain antibody fragments.](#)

For research use only. Not intended for human and animal therapeutic or diagnostic use.

生产商：南京金斯瑞生物科技有限公司 江苏省南京市江宁区科学园雍熙路28号

Manufacturer: Nanjing GenScript Biotech Co., Ltd. No. 28Yongxi Road, Jiangning District, Nanjing, Jiangsu, China