

## DATASHEET

Version: 201807

### 10 x *Taq* Buffer

**Cat. No.** B0005

**Name:** 10 x *Taq* Buffer

**Size:** 1.5 mL

**Description:** The 10X *Taq* Buffer contains 500 mM KCl, 100 mM Tris-HCl (pH 9.0 at 25°C), 15 mM MgCl<sub>2</sub>, and 1% Triton X-100. This buffer is optimized for use with 200 μM dNTPs.

**Storage:** Store at -20°C.

**Application:** Applied for PCR amplification.