

CERTIFICATE OF ANALYSIS

Product Information

Product Name:	Fms-related tyrosine kinase 3 ligand (Flt-3L), Human
Cat. No.:	Z02926-10; Z02926-50; Z02926-1
Lot No.:	B60031208
Size:	10µg/vial; 50µg /vial; 1mg/vial
Storage Condition:	Lyophilized recombinant human Fms-related tyrosine kinase 3 ligand (Flt-3L) remains stable up to 6 months at lower than -70°C from date of receipt. Upon reconstitution, hFlt-3L should be stable up to 1 week at 4°C or up to 2 months at -20°C.

Quality Control

Test Items	Specifications	Results
Physical Appearance	White Freeze-dried Powder	Conformed
Purity (SDS-PAGE)	> 95%	Conformed
Purity (HPLC)	> 95%	Conformed
Specific Activity	< 1 ng/ml	0.38 ng/ml
Molecular Weight	18~22 kDa	Conformed
Bacterial Endotoxin	< 0.2 EU/µg	Conformed

Caution:

For laboratory use only. Not for household use. If you have any questions about the Certificate of Analysis, please contact our customer service representative at 1-877-436-7274 (Toll-Free), or 1-732-885-9188.

Certified by:



Date:08/30/2012

Bioprocess Department