

DATASHEET

Version: 201807

dATP (100 mM) ≥99% HPLC**Cat. No.** C01577-10**Name:** dATP (100 mM) ≥99% HPLC**Synonyms:** 2'-Deoxyadenosine-5'-Triphosphate, Sodium Salt**Size:** 10 mL(10.0mL/vial)**Description:** dATP (100 mM) ≥99% HPLC.**Cas No.:** 1927-31-7**Molecular Formula:** $C_{10}H_{13}N_5Na_3O_{12}P_3$ **MW:** 557.2**Storage:** Stable at -20°C. Always avoid freeze-thaw cycles or exposure to frequent temperature changes.