

**dATP (100 mM) ≥99% HPLC****DATASHEET**

Version: 2016-07-19

**Cat. No.** C01577-1

cycles or exposure to frequent temperature changes.

**Name:** dATP (100 mM) ≥99% HPLC**Synonyms:** 2'-Deoxyadenosine-5'-Triphosphate,  
Sodium Salt**Size:** 1 ml(1.0ml/vial)**Description:** dATP (100 mM) ≥99% HPLC**Cas No.:** 1927-31-7**Molecular Formula:** C<sub>10</sub>H<sub>13</sub>N<sub>5</sub>Na<sub>3</sub>O<sub>12</sub>P<sub>3</sub>**MW:** 557.2**Storage:** Stable at -20°C. Always avoid freeze-thaw