

dATP (100 mM) ≥99% HPLC**DATASHEET**

Version: 2016-12-16

Cat. No. C01577-250

cycles or exposure to frequent temperature changes.

Name: dATP (100 mM) ≥99% HPLC**Synonyms:** 2'-Deoxyadenosine-5'-Triphosphate,
Sodium Salt**Size:** 250 ml (250.0ml/bottle)**Description:** dATP (100 mM) ≥99% HPLC**Cas No.:** 1927-31-7**Molecular Formula:** C₁₀H₁₃N₅Na₃O₁₂P₃**MW:** 557.2**Storage:** Stable at -20°C. Always avoid freeze-thaw