

dGTP (100 mM) ≥99% HPLC**DATASHEET**

Version: 2016-12-16

Cat. No. C01578-250

cycles or exposure to frequent temperature changes.

Name: dGTP (100 mM) ≥99% HPLC**Synonyms:** 2'-Deoxyguanosine-5'-Triphosphate,
Sodium Salt**Size:** 250 ml (250.0ml/bottle)**Description:** dGTP (100 mM) ≥99% HPLC**Cas No.:** 93919-41-6**Molecular Formula:** C₁₀H₁₃N₅Na₃O₁₃P₃**MW:** 573.1**Storage:** Stable at -20°C. Always avoid freeze-thaw