
DATASHEET

Version: 201807

dNTP (25 mM each) ≥99% HPLC**Cat. No.** C01581-250**Name:** dNTP (25 mM each) ≥99% HPLC**Size:** 250 mL(250.0mL/bottle)

Description: 25 mM dNTP Mix is suitable for use in PCR, sequencing, nick translation, fill-in, cDNA synthesis, TdT tailing reactions, and dilution of radiolabeled dNTPs. These nucleotides are supplied as convenient, ready-to-use solutions in sterilized water. The 25 mM dNTP Mix contains all four of the dNTPs, each at a final concentration of 25 mM.

Storage: Stable at -20°C. Always avoid freeze-thaw cycles or exposure to frequent temperature changes.