

## DATASHEET

Version: 201807

### **dNTP (10 mM each) $\geq$ 99% HPLC**

**Cat. No.** C01582-1

**Name:** dNTP (10 mM each)  $\geq$ 99% HPLC

**Size:** 1 mL(1.0mL/vial)

**Description:** 10 mM dNTP Mix is suitable for use in PCR, sequencing, nick translation, fill-ins, cDNA synthesis, TdT tailing reactions, and the dilution of radiolabeled dNTPs. These nucleotides are supplied as convenient, ready-to-use solutions in sterilized water. The 10 mM dNTP Mix contains all four dNTPs, each at a final concentration of 10 mM.

**Storage:** Stable at  $-20^{\circ}\text{C}$ . Always avoid freeze-thaw cycles or exposure to frequent temperature changes.