

DATASHEET

Version: 201807

dNTP (10 mM each) ≥99% HPLC**Cat. No.** C01582-250**Name:** dNTP (10 mM each) ≥99% HPLC**Size:** 250 mL(250.0mL/bottle)

Description: 10 mM dNTP Mix is suitable for use in PCR, sequencing, nick translation, fill-ins, cDNA synthesis, TdT tailing reactions, and the dilution of radiolabeled dNTPs. These nucleotides are supplied as convenient, ready-to-use solutions in sterilized water. The 10 mM dNTP Mix contains all four dNTPs, each at a final concentration of 10 mM.

Storage: Stable at -20°C. Always avoid freeze-thaw cycles or exposure to frequent temperature changes.