

Rev01  
Update: Sep,04,2025

**DATASHEET**

# dNTP (10 mM each) $\geq$ 99% HPLC

Cat. No.: C01582

## Product Introduction

<b>Product Description</b>	10 mM dNTP Mix is suitable for use in PCR, sequencing, nick translation, fill-ins, cDNA synthesis, TdT tailing reactions, and the dilution of radiolabeled dNTPs. These nucleotides are supplied as convenient, ready-to-use solutions in sterilized water. The 10 mM dNTP Mix contains all four dNTPs, each at a final concentration of 10 mM
<b>Storage</b>	Stable at -20°C. Always avoid freeze-thaw cycles or exposure to frequent temperature changes.

**For laboratory research use only. Direct human use, including taking orally and injection and clinical use are forbidden.**

Manufacturer: Nanjing GenScript Biotech Co., Ltd. No. 28Yongxi Road, Jiangning District, Nanjing, Jiangsu, China