

DATASHEET

Version: 201807

dUTP (100 mM) ≥99% HPLC**Cat. No.** C01717-1; C01717-10; C01717-50; C01717-250**Name:** dUTP (100 mM) ≥99% HPLC**Synonyms:** 2'-Deoxyadenosine-5'-Triphosphate, Sodium Salt**Size:** 1mL; 10mL; 50mL; 250mL**Description:** dUTP (100 mM) ≥99% HPLC.**Cas No.:** 102814-08-4**Molecular Formula:** C₉H₁₂N₂Na₃O₁₄P₃**MW:** 534.1**Storage:** Stable at -20°C. Always avoid freeze-thaw cycles or exposure to frequent temperature changes.