


# Mass Spectrum



## Sample Information

User : Gary  
Month-Day Processed : 03/11/16  
Time Processed : 15:08:30  
Injection Volume : 0.2  
Sample Name : Influenza Hemagglutinin (307-319)  
Cat No. : PE9051603  
Lot No. : RP20087  
Theoretical MW : 1503.79  
Observed MW : 1503.60

Interface : ESI  
Nebulizing Gas Flow : 1.5L/min  
CDL Temp : 250°C  
Block Temp : 200°C

Interface Bias : +4.5kV  
Drying Gas Flow : 5 L/min  
T.Flow : 0.2ml/min  
B.conc : 50%H2O/50%MeOH