

Sample Name : β -amyloid (1-11)
 Cat. No. : RP20106
 Lot. No. : PE3861604


Pump A : 0.065% trifluoroacetic in 100% water (v/v)
 Pump B : 0.05% trifluoroacetic in 100% acetonitrile (v/v)
 Total Flow: 1 ml/min
 Wavelength: 220 nm

Time	Unit	Command	Value	Comment
0.01	Pumps	Pump B Conc.	5	
25.00	Pumps	Pump B Conc.	65	
25.01	Pumps	Pump B Conc.	95	
31.00	Pumps	Pump B Conc.	95	
31.01	Pumps	Pump B Conc.	5	
40.00	Pumps	Pump B Conc.	5	
40.01	Controller	Stop		

<<Column Performance>>

<Detector A>

Column : AlltimaTM C18 4.6 x 250 mm



1 Det. A Ch1 / 220nm

Peak Table

Detector A Ch1 220nm

Peak#	Ret. Time	Area	Height	Area %
1	9.875	5016	1104	0.342
2	10.100	1394317	256773	95.073
3	10.352	19745	3129	1.346
4	10.554	702	198	0.048
5	10.806	23971	3971	1.635
6	10.900	1962	1013	0.134
7	12.249	20861	880	1.422
Total		1466573	267068	100.000