

Sample Name : β -Amyloid (1-42), rat
 Cat. No. : RP10013
 Lot. No. : PE1731601

Pump A : 0.065% ammonia water in 100% water (v/v)
 Pump B : 0.05% ammonia water in 100% acetonitrile (v/v)
 Total Flow: 1 ml/min
 Wavelength: 220 nm


Time	Unit	Command	Value	Comment
0.01	Pumps	Pump B Conc.	5	
25.00	Pumps	Pump B Conc.	65	
25.01	Pumps	Pump B Conc.	95	
31.00	Pumps	Pump B Conc.	95	
31.01	Pumps	Pump B Conc.	5	
37.00	Pumps	Pump B Conc.	5	
37.01	Controller	Stop		

<<Column Performance>>

<Detector A>

Column : UniPS 10-300 4.6 x 250 mm

Chromatogram



1 Det. A Ch1 / 220nm

Peak Table

Detector A Ch1 220nm

Peak#	Ret. Time	Area	Height	Area %
1	12.125	200692	11020	1.099
2	12.742	17445263	563139	95.549
3	19.376	437256	14426	2.395
4	21.984	46517	2915	0.255
5	22.847	128225	8205	0.702
Total		18257953	599705	100.000