


Mass Spectrum



Sample Information

User : Gary
Month-Day Processed : 03/29/16
Time Processed : 16:28:22
Injection Volume : 0.2
Sample Name : Influenza Hemagglutinin (307-319)
Cat No. : PE0841603
Lot No. : RP20087
Theoretical MW : 1503.79
Observed MW : 1503.45

Interface : ESI
Nebulizing Gas Flow : 1.5L/min
CDL Temp : 250°C
Block Temp : 200°C

Interface Bias : +4.5kV
Drying Gas Flow : 5 L/min
T.Flow : 0.2ml/min
B.conc : 50%H2O/50%MeOH