


# Mass Spectrum



**Sample Information**  
User : Gary  
Month-Day Processed : 07/13/16  
Time Processed : 14:04:36  
Injection Volume : 0.2  
Sample Name : Gastrin-1, human  
Cat No. : RPI2740  
Lot No. : PE7601606  
Theoretical MW : 2098.21  
Observed MW : 2097.20

**Interface**  
Interface : ESI  
Nebulizing Gas Flow : 1.5L/min  
CDL Temp : 250°C  
Block Temp : 200°C

**Interface Bias**  
Drying Gas Flow : +4.5kV  
T.Flow : 5 L/min  
B.conc : 0.2ml/min  
:50%H2O/50%MeOH