


Mass Spectrum



Sample Information

User : Gary
Month-Day Processed : 07/26/16
Time Processed : 10:47:16
Injection Volume : 0.2
Sample Name : LCMV_gp33-41
Cat No. : RP20257
Lot No. : PE1711607
Theoretical MW : 1044.25
Observed MW : 1043.40

Interface : ESI
Nebulizing Gas Flow : 1.5L/min
CDL Temp : 250°C
Block Temp : 200°C

Interface Bias : +4.5kV
Drying Gas Flow : 5 L/min
T.Flow : 0.2ml/min
B.conc : 50%H₂O/50%MeOH