

Sample Name : BDC2.5 mimotope 1040-31  
 Cat. No. : RP20518  
 Lot. No. : PE9801609

Pump A : 0.065% trifluoroacetic in 100% water (v/v)  
 Pump B : 0.05% trifluoroacetic in 100% acetonitrile (v/v)  
 Total Flow:1 ml/min  
 Wavelength:220 nm


Time	Unit	Command	Value	Comment
0.01	Pumps	Pump A B.Conc	5	
25.00	Pumps	Pump A B.Conc	65	
25.01	Pumps	Pump A B.Conc	95	
31.00	Pumps	Pump A B.Conc	95	
31.01	Pumps	Pump A B.Conc	5	
40.00	Pumps	Pump A B.Conc	5	
40.01	Controller	Stop		

<<Column Performance>>

<Detector A>

Column : Alltima™ C18 4.6 x 250 mm

Chromatogram



1 Det.A Ch1 / 220nm

Peak Table

Detector A Ch1 220nm

Peak#	Ret. Time	Area	Height	Area %
1	16.092	41132	5862	0.227
2	17.390	17945269	1898543	99.245
3	17.741	70775	10955	0.391
4	17.963	7370	1151	0.041
5	18.066	8477	1395	0.047
6	18.923	8845	1531	0.049
Total		18081867	1919436	100.000