

DATASHEET Version: 201807

dNTP mixture, 10 mM each

Cat. No. D0056

Name: dNTP mixture, 10 mM each

Size: 0.5 mL(0.5 mL/vial)

Description: 10 mM dNTP Mix is suitable for use in PCR, sequencing, nick translation, fill-ins, cDNA synthesis, TdT tailing reactions, and the dilution of radiolabeled dNTPs. These nucleotides are supplied as convenient, ready-to-use solutions in sterilized water. The 10 mM dNTP Mix contains all four dNTPs, each at a final concentration of 10 mM.

Storage: Stable at -20°C. Always avoid freeze-thaw cycles or exposure to frequent temperature changes.