

IFN-tau, Ovine

Cat. No.: Z02790-10

Size: 10.0 ug

Synonyms: Interferon-tau (IFN-tau), Ovine;

Description:

IFN-tau is a new class of type I IFN that is secreted by the trophoblast and is the signal for maternal recognition of pregnancy in sheep. IFN- tau has potent immunosuppressive and antiviral activities similar to other type I IFN but is less cytotoxic than IFN- α/β . The current investigation concerns the effect of recombinant ovine IFN- tau (rOvIFN- tau) on the modulation of MHC class I and II expression on cloned mouse cerebrovascular endothelial (CVE) cells. IFN- tau induced tyrosine phosphorylation of Stat1 and up-regulated the expression of MHC class I on CVE. One proposed action by which type I IFN reduce the relapse rate in MS is via interference with IFN- γ -induced MHC class II expression. IFN- tau was shown to downregulate IFN- γ -induced MHC class II expression on CVE and, hence, may be of potential therapeutic value in downregulating inflammation in the central nervous system (CNS). IFN- tau did not upregulate the expression of MHC class II on CVE. IFN- tau also inhibited the replication of Theiler's virus in CVE.

Amino Acid Sequence:

```
00001 CYLSRKLMLD ARENLKLLDR MNRLSPHSCL QDRKDFGLPQ
00041 EMVEGDQLQK DQAFPVLYEM LQQSFNLFYT EHSSAAWDTT
00081 LLEQLCTGLQ QQLDHLDTCR GQVMGEEDSE LGNMDPIVTV
00121 KKYFQGIYDY LQEKGYSDCA WEIVRVEMMR ALTVSTTLQK
00161 RLTKMGDDLN SP
```

Source: *P. pastoris*

Species: Ovine

Biological Activity: Fully biologically active when compared to IFN-alpha. The specific activity determined by a viral resistance assay is no less than 1.0×10^7 IU/mg.

Molecular Weight: Approximately 19.9 kDa, a single non-glycosylated polypeptide chain containing 172 amino acids.

Formulation: Lyophilized from a 0.2 μ m filtered concentrated solution in PBS, pH 7.4.

Appearance: Sterile Filtered White lyophilized (freeze-dried) powder.

Reconstitution: We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Reconstitute in sterile distilled water or aqueous buffer containing 0.1 % BSA to a concentration of 0.1-1.0 mg/mL. Stock solutions should be apportioned into working aliquots and stored at $\leq -20^\circ\text{C}$. Further dilutions should be made in appropriate buffered solutions.

Purity: > 95 % by SDS-PAGE and HPLC analyses.

Endotoxin Level: Less than 1 EU/ μ g of rOvIFN- τ as determined by LAL method.

Storage: This lyophilized preparation is stable at 2-8 $^\circ\text{C}$, but should be kept at -20°C for long term storage, preferably desiccated. Upon reconstitution, the preparation is stable for up to one week at 2-8 $^\circ\text{C}$. For maximal stability, apportion the reconstituted preparation into working aliquots and store at -20°C to -70°C . Avoid repeated freeze/thaw cycles.