

HBsAg, adr

Cat. No.: Z02813-1

Size: 1.0 mg

Synonyms: Hepatitis B Surface Antigen adr Subtype (HBsAg, adr)

Description:

Hepatitis B Surface Antigen adr Subtype(HBsAg, adr)

Source: CHO

Biological Activity: Was tested in ELISA with anti-HBsAg antibodies.

Molecular Weight: Approximately 23.0-27.0 kDa on a reducing SDS-PAGE.

Formulation: Sterile Filtered solution containing 20 mM PB, 154 mM sodium chloride.

Appearance: Sterile Filtered colourless liquid

Purity: >95% by HPLC analyses.

Endotoxin Level: Less than 0.2EU/ug of rHBsAgadr as determined by LAL method.

Storage: Recombinant Hepatitis B Surface Antigen can be stored under sterile conditions at 2-8°C for up to six months. DO NOT FREEZE.