

HBsAg, adr

Cat. No.: Z02813-100

Size: 100.0 ug

Synonyms: Hepatitis B Surface Antigen adr Subtype (HBsAg, adr)

Description:

Hepatitis B Surface Antigen adr Subtype(HBsAg, adr)

Source: CHO

Biological Activity: Was tested in ELISA with anti-HBsAg antibodies.

Molecular Weight: Approximately 23.0-27.0 kDa on a reducing SDS-PAGE.

Formulation: Sterile Filtered solution containing 20 mM PB, 154 mM sodium chloride.

Appearance: Sterile liquid.

Reconstitution: We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom.

Purity: >95% by HPLC analyses.

Endotoxin Level: Less than 0.2EU/ug of rHBsAgadr as determined by LAL method.

Storage: Recombinant Hepatitis B Surface Antigen can be stored under sterile conditions at 2-8°C for up to six months. DO NOT FREEZE.