

HBsAg, adw

Cat. No.: Z02814-1

Size: 1.0 mg

Synonyms: Hepatitis B Surface Antigen adw Subtype (HBsAg, adw)

Description:

Hepatitis B Surface Antigen adw Subtype (HBsAg, adw)

Source: *P. pastoris*

Biological Activity: Was tested in ELISA with anti-HBsAg antibodies

Molecular Weight: Approximately 24.0 kDa on a reducing SDS-PAGE.

Formulation: Sterile Filtered solution containing 20 mM PB, 154 mM sodium chloride.

Appearance: Sterile Filtered colourless liquid

Purity: >95% by reducing SDS-PAGE and HPLC analyses

Endotoxin Level: Less than 0.2EU/ug of rHBsAgadw as determined by LAL method.

Storage: Recombinant Hepatitis B Surface Antigen can be stored under sterile conditions at 2-8°C for up to six months. DO NOT FREEZE.