

## HBsAg, adw

**Cat. No.:** Z02814-100

**Size:** 100.0 ug

**Synonyms:** Hepatitis B Surface Antigen adw Subtype (HBsAg, adw)

**Description:**

Hepatitis B Surface Antigen adw Subtype (HBsAg, adw)

**Source:** *P. pastoris*

**Biological Activity:** Was tested in ELISA with anti-HBsAg antibodies

**Molecular Weight:** Approximately 24.0 kDa on a reducing SDS-PAGE.

**Formulation:** Sterile Filtered solution containing 20 mM PB, 154 mM sodium chloride.

**Appearance:** Sterile Filtered colourless liquid

**Purity:** >95% by reducing SDS-PAGE and HPLC analyses

**Endotoxin Level:** Less than 0.2EU/ug of rHBsAgadw as determined by LAL method.

**Storage:** Recombinant Hepatitis B Surface Antigen can be stored under sterile conditions at 2-8°C for up to six months. DO NOT FREEZE.