

Sample Name :KRAS G12C Peptide (VVVGACGVGK)
Cat No. : RP30236
Lot No. : E0033626I
Time Processed :20:56:32
Month-Day-Year Processed :09/03/2026

Pump A : 0.065% trifluoroacetic in 100% water (v/v)
Pump B : 0.05% trifluoroacetic in 100% acetonitrile (v/v)
Total Flow:1 ml/min
Wavelength:220 nm

<<LC Time Program>>

Time	Module	Command	Value
0.01	Pumps	B.Conc	5
25.00	Pumps	B.Conc	65
25.01	Pumps	B.Conc	95
27.00	Pumps	B.Conc	95
27.01	Pumps	B.Conc	5
35.00	Pumps	B.Conc	5
35.01	Controller	Stop	

<<Column Performance>>


<Detector A>

Column :Inertsil ODS-3 4.6 x 250 mm

Equipment: ZJ19010325

<Chromatogram>

mV



<Peak Table>

Detector A 220nm

Peak#	Ret. Time	Area	Height	Area%
1	8.350	8806	993	0.351
2	8.717	33616	5231	1.342
3	8.912	2420082	332326	96.588
4	9.252	17850	2181	0.712
5	9.711	16096	2693	0.642
6	10.798	9123	558	0.364
Total		2505573	343982	100.000