

Sample Name :Melittin,honey bee
Cat No. : RP20415
Lot No. : E0124624J
Time Processed :18:01:25
Month-Day-Year Processed :10/10/2024

Pump A : 0.065% trifluoroacetic in 100% water (v/v)
Pump B : 0.05% trifluoroacetic in 100% acetonitrile (v/v)
Total Flow:1 ml/min
Wavelength:220 nm

<<LC Time Program>>

Time	Module	Command	Value
0.01	Pumps	B.Conc	5
25.00	Pumps	B.Conc	65
25.01	Pumps	B.Conc	95
27.00	Pumps	B.Conc	95
27.01	Pumps	B.Conc	5
35.00	Pumps	B.Conc	5
35.01	Controller	Stop	

<<Column Performance>>


<Detector A>

Column :Inertsil ODS-SP 4.6 x 250 mm

Equipment: ZJ20010140

<Chromatogram>

mV



<Peak Table>

Detector A 220nm

Peak#	Ret. Time	Area	Height	Area%
1	9.583	20690	2474	0.169
2	11.142	18021	2551	0.147
3	12.942	11748362	1018527	96.013
4	13.367	402829	44746	3.292
5	13.567	35426	5005	0.290
6	13.700	10874	2436	0.089
Total		12236203	1075739	100.000