

Sample Name :Melittin,honey bee  
 Cat No. : RP20415  
 Lot No. : E0155422J  
 Time Processed :21:15:36  
 Month-Day-Year Processed :11/01/2022  
 Pump A : 0.065% trifluoroacetic in 100% water (v/v)  
 Pump B : 0.05% trifluoroacetic in 100% acetonitrile (v/v)  
 Total Flow:1 ml/min  
 Wavelength:220 nm

<<LC Time Program>>

| Time  | Module     | Command | Value |
|-------|------------|---------|-------|
| 0.01  | Pumps      | B.Conc  | 5     |
| 25.00 | Pumps      | B.Conc  | 65    |
| 25.01 | Pumps      | B.Conc  | 95    |
| 27.00 | Pumps      | B.Conc  | 95    |
| 27.01 | Pumps      | B.Conc  | 5     |
| 35.00 | Pumps      | B.Conc  | 5     |
| 35.01 | Controller | Stop    |       |

<<Column Performance>>


<Detector A>

Column :Inertsil ODS-SP 4.6 x 250 mm

Equipment: ZJ19010325

### <Chromatogram>

mV



### <Peak Table>

Detector A Channel 1 220nm

| Peak# | Ret. Time | Area     | Height  | Area%   |
|-------|-----------|----------|---------|---------|
| 1     | 19.609    | 29060    | 2685    | 0.189   |
| 2     | 20.043    | 411807   | 20969   | 2.685   |
| 3     | 20.414    | 14603756 | 1218847 | 95.204  |
| 4     | 20.746    | 266387   | 27477   | 1.737   |
| 5     | 21.420    | 28442    | 2639    | 0.185   |
| Total |           | 15339452 | 1272617 | 100.000 |