


Mass Spectrum



Sample Information

Month-Day Processed : 08/09/23

Time Processed : 12:45:24

Injection Volume : 0.2

Sample Name : C-Type Natriuretic Peptide (CNP) (1-22), human

Cat No. : RP11110

Lot No. : E0170523G

Theoretical MW : 2197.60

Observed MW : 2197.8

Interface

Nebulizing Gas Flow : 1.5L/min

CDL Temp : 250

Block Temp : 200

Equipment

Interface Bias : +4.5 kV

Drying Gas Flow : 5 L/min

T.Flow : 0.2 ml/min

B.conc : 50%H₂O/50%MeOH