


Mass Spectrum



Sample Information

Month-Day Processed :	02/03/26	Equipment	: ZJ21010035
Time Processed :	23:20:49	Interface Bias	: +4.5 kV
Injection Volume :	0.8	Drying Gas Flow	: 5 L/min
Sample Name :	瑞淑颜™类蛇毒肽	T.Flow	: 0.2 ml/min
Cat No. :	RP30094CN	B.conc	: 50%H ₂ O/50%MeOH
Lot No. :	E0175026A		
Theoretical MW :	495.6		
Observed MW :	375.1		