

Sample Name :Hepatitis B Virus Core (128-140)
 Cat No. :RP17313
 Lot No. :E0203422H
 Time Processed :6:15:23 PM
 Month-Day-Year Processed :08/15/2022

Pump A : 0.065% trifluoroacetic in 100% water (v/v)
 Pump B : 0.05% trifluoroacetic in 100% acetonitrile (v/v)
 Total Flow:1 ml/min
 Wavelength:220 nm

<<LC Time Program>>

Time	Module	Command	Value
0.01	Pumps	B.Conc	5
25.00	Pumps	B.Conc	65
25.01	Pumps	B.Conc	95
27.00	Pumps	B.Conc	95
27.01	Pumps	B.Conc	5
35.00	Pumps	B.Conc	5
35.01	Controller	Stop	


<<Column Performance>>

<Detector A>

Column :Inertsil ODS-3 4.6 x 250 mm
 Equipment: ZJ19010141

<Chromatogram>

mV



<Peak Table>

Detector A Channel 1 220nm

Peak#	Ret. Time	Area	Height	Area%
1	10.722	23903	2544	0.094
2	11.642	11405	901	0.045
3	12.318	385793	29750	1.515
4	12.787	24903299	1819948	97.780
5	13.497	91165	5623	0.358
6	14.986	53209	2768	0.209
Total		25468774	1861534	100.000