


Mass Spectrum



Sample Information

Month-Day Processed :	11/09/25	Equipment :	ZJ21010035
Time Processed :	13:01:47	Interface :	ESI
Injection Volume :	0.3	Nebulizing Gas Flow :	1.5L/min
Sample Name :	C-Type Natriuretic Peptide (CNP) (1-22), human	CDL Temp :	250
Cat No. :	RP11110	Block Temp :	200
Lot No. :	E0353625K		
Theoretical MW :	2197.60		
Observed MW :	2197.2		

Interface :	ESI	Equipment :	ZJ21010035
Nebulizing Gas Flow :	1.5L/min	Interface Bias :	+4.5 kV
CDL Temp :	250	Drying Gas Flow :	5 L/min
Block Temp :	200	T.Flow :	0.2 ml/min
		B.conc :	50%H ₂ O/50%MeOH