


Mass Spectrum



Sample Information		Equipment : ZI21010035	
Month-Day Processed :	06/20/23	Interface Bias :	+4.5 kV
Time Processed :	19:41:45	Drying Gas Flow :	5 L/min
Injection Volume :	0.3	T.Flow :	0.2 ml/min
Sample Name :	DYKDDDDK Peptide	B. conc :	50% H2O/50% MeOH
Cat No. :	RP10586		
Lot No. :	E0491123F		
Theoretical MW :	1012.97		
Observed MW :	1013.0		