


Mass Spectrum



Sample Information

Month-Day Processed :	06/27/23	Equipment	:ZJ22010150
Time Processed :	20:16:19	Interface Bias	: +4.5 kV
Injection Volume :	0.2	Drying Gas Flow	:5 L/min
Sample Name :	TAT peptide	T.Flow	:0.2 ml/min
Cat No. :	RP20256	B.conc	:50%H2O/50%MeOH
Lot No. :	E0514623F		
Theoretical MW :	1621.91		
Observed MW :	1622.0		