

Sample Name :C-Type Natriuretic Peptide (CNP) (1-22), human
 Cat No. :RP11110
 Lot No. :E0519122G
 Time Processed : 17:53:29
 Month-Day-Year Processed : 07/27/2022

Pump A : 0.065% trifluoroacetic in 100% water (v/v)
 Pump B : 0.05% trifluoroacetic in 100% acetonitrile (v/v)
 Total Flow:1 ml/min
 Wavelength:220 nm

Time	Module	Command	Value
0.01	Pumps	Pump A B.Conc	5
25.00	Pumps	Pump A B.Conc	65
25.01	Pumps	Pump A B.Conc	95
27.00	Pumps	Pump A B.Conc	95
27.01	Pumps	Pump A B.Conc	5
35.00	Pumps	Pump A B.Conc	5
35.01	Controller	Stop	


<<Column Performance>>

<Detector A>

Column : Inertsil ODS-3 4.6 x 250 mm

Equipment: GK11010017

mV



1 Detector A Channel 1 / 220nm

Peak Table

Detector A Channel 1 220nm

Peak#	Ret. Time	Area	Height	Area %
1	10.606	8051	1277	0.071
2	13.670	12545	1535	0.110
3	13.859	8461	902	0.074
4	14.136	14382	1335	0.126
5	14.722	179870	18530	1.576
6	15.174	11040112	1153022	96.755
7	15.809	43251	4027	0.379
8	17.092	103686	8781	0.909
Total		11410358	1189409	100.000