


Mass Spectrum



Sample Information

|                          |                                                |                     |           |                 |              |
|--------------------------|------------------------------------------------|---------------------|-----------|-----------------|--------------|
| Month-Day Processed :    | 07/26/22                                       | Interface           | :ESI      | Equipment       | : ZJ21010035 |
| Time Processed :         | 15:10:31                                       | Nebulizing Gas Flow | :1.5L/min | Interface Bias  | : +4.5 kV    |
| Injection Volume :       | 0.4                                            | CDL Temp            | :250      | Drying Gas Flow | :5 L/min     |
| Sample Name :            | C-Type Natriuretic Peptide (CNP) (1-22), human | Block Temp          | :200      | T.Flow          | :0.2 ml/min  |
| Cat No. :                | RP11110                                        |                     |           | B.conc          | :50%H2O/50%M |
| Lot No. :                | E0519122G                                      |                     |           |                 |              |
| Theoretical MW : 2197.60 |                                                |                     |           |                 |              |
| Observed MW : 2197.2     |                                                |                     |           |                 |              |