


Mass Spectrum



Sample Information

Month-Day Processed :	02/17/22	Equipment	: GK11010007
Time Processed :	1:23:24 PM	Interface Bias	: +4.5 kV
Injection Volume :	0.4	Drying Gas Flow	: 5 L/min
Sample Name :	Amylin, human, amide	T.Flow	: 0.2 ml/min
Cat No. :	RP11278	B.conc	: 50% H ₂ O/50% MeOH
Lot No. :	E0591422B		
Theoretical MW :	3903.28		
Observed MW :	3903.2		