


Mass Spectrum



Sample Information

|                         |                   |                 |                 |
|-------------------------|-------------------|-----------------|-----------------|
| Month-Day Processed :   | 10/13/24          | Equipment       | :ZJ22010150     |
| Time Processed :        | 18:39:03          | Interface Bias  | : +4.5 kV       |
| Injection Volume :      | 0.3               | Drying Gas Flow | :5 L/min        |
| Sample Name :           | Melitin,honey bee | T.Flow          | :0.2 ml/min     |
| Cat No. : RP20415       |                   | B.conc          | :50%H2O/50%MeOH |
| Lot No. : E0596524J     |                   |                 |                 |
| Theoretical MW: 2847.50 |                   |                 |                 |
| Observed MW : 2847.3    |                   |                 |                 |