

Sample Name :OVA Peptide (323-339)
Cat No. :RP10610
Lot No. :E0637725L
Time Processed :14:15:41
Month-Day-Year Processed :12/23/2025
Pump A : 0.065% trifluoroacetic in 100% water (v/v)
Pump B : 0.05% trifluoroacetic in 100% acetonitrile (v/v)
Total Flow:1 ml/min
Wavelength:220 nm

<<LC Time Program>>

Time	Module	Command	Value
0.01	Pumps	B.Conc	5
25.00	Pumps	B.Conc	65
25.01	Pumps	B.Conc	95
27.00	Pumps	B.Conc	95
27.01	Pumps	B.Conc	5
35.00	Pumps	B.Conc	5
35.01	Controller	Stop	

<<Column Performance>>


<Detector A>

Column :Inertsil ODS-SP 4.6 x 250 mm

Equipment: ZJ21010376

<Chromatogram>

mV



<Peak Table>

Detector A 220nm

Peak#	Ret. Time	Area	Height	Area%
1	7.817	36361	5986	0.842
2	8.208	4120059	651785	95.407
3	8.367	70585	21450	1.635
4	8.417	76580	17611	1.773
5	8.767	10737	1537	0.249
6	8.883	4063	890	0.094
Total		4318385	699261	100.000