

Sample Name :Influenza A NP(366-374) Strain A/PR/8/35

Cat No. : RP20267

Lot No. : E0669024B

Time Processed :10:43:40

Month-Day-Year Processed :02/21/2024

Pump A : 0.065% trifluoroacetic in 100% water (v/v)

Pump B : 0.05% trifluoroacetic in 100% acetonitrile (v/v)

Total Flow:1 ml/min

Wavelength:220 nm

<<LC Time Program>>

Time	Module	Command	Value
0.01	Pumps	B.Conc	5
25.00	Pumps	B.Conc	65
25.01	Pumps	B.Conc	95
27.00	Pumps	B.Conc	95
27.01	Pumps	B.Conc	5
35.00	Pumps	B.Conc	5
35.01	Controller	Stop	

<<Column Performance>>


<Detector A>

Column :Inertsil ODS-SP 4.6 x 250 mm

Equipment: GK12010012

<Chromatogram>

mV



<Peak Table>

Detector A Channel 1 220nm

Peak#	Ret. Time	Area	Height	Area%
1	5.200	15581	2240	0.807
2	5.725	2060	292	0.107
3	6.208	1897792	182556	98.312
4	7.333	14953	1742	0.775
Total		1930386	186829	100.000